



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Synchronous Machine Refurbishing

**Product Code :**  
LBNY-0005-121000184



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Refurbishing Work  
Console Mechanical Reassembling  
Console reassembling  
Clear rear panel changing  
Assembling and mechanical adjustment of motor  
Repair and/or replacement of casters and couplings  
Working table replacement

Complete Console Disassembling  
Mechanical repair of housing if required  
Sand blasting and painting

Complete Testing And Verifying  
Coupling verifying on both sides of machine  
Electrical testing w/o load  
Inspection report and quality control.

Complete Motor Disassembling  
Motor inspection  
Repair and replacement of brush assembly if required  
External painting  
Bearing replacement

---

Electrical Reassembling And Rewiring  
Electrical part replacement  
Front panel replacement  
Rewiring.

**Technical Specification :**

Synchronous Machine Refurbishing



**Laboratory Instrument India**