



Email : sales@laboratoryinstrumentindia.com

Product Name :

Diesel Engine Bench Detroit Diesel 14.7L

Product Code :

LBNY-0005-121000192



www.laboratoryinstrumentindia.com

Description :

Especially built to meet all diesel engine performance and repair tasks. Completely operational and safe, this standard trainer is easy to work on and complies with manufacturer rules and specifications related to all OEM codes, colors and electronics. The Detroit Diesel engine bench w/EGR, DPF, SCR, is mounted on a tubing steel stand.

Educational Advantages:

Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation.

Demonstration of major engine systems to groups of students without the access limitations of a complete vehicle.

Technical Specification :

Diesel Engine Bench Detroit Diesel 14.7L



Laboratory Instrument India

