



Email : sales@laboratoryinstrumentindia.com

Product Name :
6.6l Diesel Engine Bench

Product Code :
LBNY-0005-121000198



Description :

The trainer is easy to work on and all major components, wiring and connectors match OEM configuration. Trainer is mounted on a mobile stand constructed of 2 (5 cm) tubular steel equipped with four 5 (12.5 cm) phenol casters. The 6.6l Diesel Engine Bench is a fully operational engine mounted on a mobile stand. Trainer is equipped with safety guards covering all rotating components. Students can perform many related Engine Performance and Electrical/Electronic tasks using this trainer.

Educational Advantages:

Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation.

Allows demonstration of engine systems to groups of students without the access limitations of a complete vehicle.

Equipped with all OEM equipment related to Emission Standards.

Technical Specification :

6.6l Diesel Engine Bench



Laboratory Instrument India

