



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Diesel Engine Bench Cummins TIER 4I

**Product Code :**  
LBNY-0005-121000201



**Description :**

The Cummins Diesel Engine Bench is mounted on a tubing steel stand. Especially built to meet all diesel engine performance and repair tasks. Completely operational and safe, this standard trainer is easy to work on and complies with manufacturer rules and specifications related to all OEM codes, colors and electronics.

**Educational Advantages:**

Engine systems respond to inserted faults with real world symptoms, OEM DTC and check engine light operation.

Demonstration of major engine systems to groups of students without the access limitations of a complete vehicle.

**Technical Specification :**

Diesel Engine Bench Cummins TIER 4I



**Laboratory Instrument India**