



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dc Machine Refurbishing

Product Code :
LBNY-0005-121000202



www.laboratoryinstrumentindia.com

Description :

Refurbishing Work
Electrical Reassembling And Rewiring
Front panel replacement
Rewiring
Electrical part replacement

Complete Testing And Verifying
Electrical testing w/o load
Coupling verifying on both sides of machine
Inspection report and quality control

Complete Console Disassembling
Mechanical repair of housing if required
Sand blasting and painting

Console Mechanical Reassembling
Console reassembling
Repair and/or replacement of casters and couplings
Working table replacement
Clear rear panel changing
Assembling and mechanical adjustment of motor

Complete Motor Disassembling

Motor inspection
Repair and replacement of brush assembly if required
Bearing replacement
External painting

Technical Specification :

Dc Machine Refurbishing



Laboratory Instrument India