



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle Lighting And Accessory Systems Trainer

Product Code :
LBNY-0005-121000208



Description :

The trainer is built to help the student evolve from building basic circuits to more complex electrical troubleshooting. Due to the fact that the trainer is fully modular, additional can be added to follow technical evolutions in motorcycle and light vehicle industry. The Motorcycle Lighting And Accessory Systems Trainer contains OEM system components (application varies) mounted on modular plastic plates called where the student completes all of the system wiring. A built-in regulated and protected 12 power supply and fault box is included with the trainer for each side. The trainer is mounted single- or double-sided to allow more than one group of students to work simultaneously on the trainer.

Educational Advantages:

The trainer can provide multiple student work stations accessible from either side
Lighting and accessory systems respond to inserted faults (12) with real world symptoms
Students can build a complete motorcycle lighting or accessories circuit and simulate a switched-off immobilizer.

Technical Specification :

Motorcycle Lighting And Accessory Systems Trainer



Laboratory Instrument India

