



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Eaton Transmission Trainer

Product Code :
LBNY-0005-121000210



Description :

This model provides several forward speeds and reverses, consisting of a front section and an auxiliary section. The auxiliary section contains LO and HI range ratios, and the L and H splitter gear ratios. The Cutaway Eaton Transmission Trainer has been designed for teaching Eaton heavy-duty transmission theory and operation. When in HI range, the gear ratios can be split by using the splitter control button. This model is driven by hand crank.

Educational Advantages:

Allows visualization of gear changes and clutch operation
Gear changes can be made using regulated air pressure and provided controls
Driven by hand crank.

Technical Specification :

Cutaway Eaton Transmission Trainer



Laboratory Instrument India

