



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Diesel Engine Bench Cummins Recycled

**Product Code :**  
LBNY-0005-121000216



**Description :**

The Recycled Cummins Diesel Engine Bench w/EGR, DPF, SCR, is mounted on a tubing steel stand. Especially built to meet all diesel engine performance and repair tasks. Completely operational and safe, this standard trainer is easy to work on and complies with manufacturer rules and specifications related to all OEM codes, colors and electronics.

**Educational Advantages:**

Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation

Functional DLC using protocols

Demonstration of major engine systems to groups of students without the access limitations of a complete vehicle.

**Technical Specification :**

Diesel Engine Bench Cummins Recycled



**Laboratory Instrument India**

