



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fuel Injection And GM Engine Management Trainer

Product Code :
LBNY-0005-121000219



Description :

The trainer is used to teach and demonstrate engine management systems in a vehicle. The Fuel Injection And GM Engine Management Trainer functions as in a real automobile, but uses non-flammable fluid instead of gasoline for safety and educational purposes.

Educational Advantages:

- Allows demonstration of engine management systems in a classroom environment
- Live ignition system function with observable spark
- Live fuel injection operation using non-flammable fluid in a transparent enclosed system
- Allows full communication with a scan tool for fault codes and data stream.

Standard Equipment and Features:

- Most engine input sensors and output actuators
- Return type fuel system with pressure regulator
- COP (Coil-On-Plug) ignition system
- All OEM components, wiring, connectors and wire colors
- OBDII data link connector
- Integrated ECM breakout box
- Adjustable sensor inputs
- ECM breakout box with test receptacles for test equipment hookup.

Technical Specification :

Fuel Injection And GM Engine Management Trainer



Laboratory Instrument India