



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop CNC Lathe

Product Code :
LBNY-0004-730001



Description :

Workshop CNC Lathe, Workshop CNC Lathe Equipments, Workshop CNC Lathe Tools, Workshop CNC Lathe Tool Kits, Workshop CNC Lathe Manufacturers, Workshop CNC Lathe Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

CNC lathe
Max Swing: 660mm
Spindle Bore: 105mm
Spindle Speed Steps Stepless in Each Step (4 Steps)
Tool Holder: 32X32mm
Main Power: 11 Kw
Application: metal
Process usage metal cutting CNC
Movement method: linear control.



Laboratory Instrument India