



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop CNC Lathe

**Product Code :**  
LBNY-0004-730001



**Description :**

Workshop CNC Lathe, Workshop CNC Lathe Equipments, Workshop CNC Lathe Tools, Workshop CNC Lathe Tool Kits, Workshop CNC Lathe Manufacturers, Workshop CNC Lathe Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

**Technical Specification :**

CNC lathe  
Max Swing: 660mm  
Spindle Bore: 105mm  
Spindle Speed Steps Stepless in Each Step (4 Steps)  
Tool Holder: 32X32mm  
Main Power: 11 Kw  
Application: metal  
Process usage metal cutting CNC  
Movement method: linear control.



**Laboratory Instrument India**