



Email : sales@laboratoryinstrumentindia.com

Product Name :
Diesel Engine Bench Cummins 2010

Product Code :
LBNY-0005-121000222



Description :

Especially built to meet all diesel engine performance and repair tasks. Completely operational and safe, this standard trainer is easy to work on and complies with manufacturer rules and specifications related to all OEM codes, colors and electronics. The Cummins 2010 Diesel Engine Bench w/EGR, DPF, SCR, is mounted on a tubing steel stand.

Educational Advantages:

Engine systems respond to inserted faults with real world symptoms, OEM DTC and check engine light operation
Demonstration of major engine systems to groups of students without the access limitations of a complete vehicle.

Technical Specification :

Diesel Engine Bench Cummins 2010



Laboratory Instrument India

