



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Power Door Lock System Trainer Four Door

**Product Code :**  
LBNY-0005-121000225



**Description :**

Terminal block receptacles are provided for diagnostic test points. The Four-Door Power Door Lock System Trainer allows for the demonstration of a complete power window system using real world components. The trainer is mounted on a plastic training board with the following features:

Electrical test points are included on the trainer for circuit diagnosis.

User manual included.

Must be connected to a fully charged automobile battery or a power supply.

OEM components (BCM, door latches, lock/unlock switches, ajar switches, etc.).

A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.

**Technical Specification :**

**Educational Advantages:**

Four operational faults can be inserted to promote student diagnostic troubleshooting competencies.

Ability to visualize the operation of power door lock latches and switches.

Ability to test and verify a basic 4-door power door lock circuit.

Includes wiring schematics.



**Laboratory Instrument India**

