



Email : sales@laboratoryinstrumentindia.com

Product Name :
Modular Electrical System Trainer

Product Code :
LBNY-0005-121000227



Description :

Assists in the teaching of Ohms Law applications as well as common electrical circuits used in automobiles
Self-contained with a built-in adjustable and protected power supply
Modules available for Ohms Law, series, parallel, combination, engine performance sensors and basic electronic circuits
Trainer is fully customizable by choosing from available modules to cover desired subject content and scope
Components available for the trainer use industry-standard schematic symbols

Main Features:

Available adjustable power supply
Protected with removable metal cover.

Technical Specification :

Modular Electrical System Trainer



Laboratory Instrument India

