



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Lathe Machine

Product Code :
LBNY-0004-730002



Description :

Workshop Lathe Machine

Technical Specification :

Lathe Machine
Machining Capacity: Medium Duty
Turned Length (mm): 1350 mm
Spindle Motor Power(kW): 5.5
Automatic Grade: Manual
Swing Diameter: 360mm over bed
Power(W): 5.5 KW
Spindle Bore(mm): 80 mm
Width of bed (mm): 336
Type: Horizontal
Dimension(L*W*H): 2810x1100x1255mm
Weight (KG): 2550 KG
CNC or Not: Normal
Length Workpiece: 1500 mm.

Workshop Lathe Machine, Workshop Lathe Machine Equipments, Workshop Lathe Machine Tools, Workshop Lathe Machine Tool Kits, Workshop Lathe Machine Manufacturers, Workshop Lathe Machine Suppliers from India, China, Kenya



Laboratory Instrument India