



Email : sales@laboratoryinstrumentindia.com

Product Name :
Can Bus Vehicle Lighting And Accessory System Trainer

Product Code :
LBNY-0005-121000230



Description :

A fully functional vehicle lighting and accessory system can be used for demonstrations and diagnosis/troubleshooting assignments.
Multi-function switch operates as in a real vehicle.
Scan tools can be connected to the DLC providing OEM data stream information, trouble codes and diagnostic PIDs for student diagnosis and troubleshooting.
Students have clear access to all vehicle lighting and accessory systems and components without the access limitations of body panels.
Operational faults can be inserted into the system for real-world experience with system problems.
Functional instrument panel with indicators
Students can use included wiring diagrams to assist in trainer diagnosis and troubleshooting
Trainer comes equipped with a 120VAC to 12VDC protected power supply
DLC allow OEM bi-directional control using a scan tool for applicable circuits
Constructed from a real vehicle, any needed repair parts can be obtained from local auto parts distributors.

Technical Specification :

Can Bus Vehicle Lighting And Accessory System Trainer

Laboratory Instrument India

