



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Custom Gasoline Direct Injection Engine Bench

**Product Code :**  
LBNY-0005-121000231



**Description :**

This custom trainer can be designed to incorporate any vehicle system you require. A Custom Gasoline Engine Bench equipped with an engine from any vehicle of your choice.

**Standard Equipment and Features:**

- All OEM components, wiring, connectors and wire colors
- Integrated ECM breakout box
- Emergency stop button
- Manual fault box for engine performance
- No component simulation
- OBDII data link connector.

**Educational Advantages:**

- Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation
- Demonstration of major vehicle systems to groups of students without the access limitations of a complete vehicle

**Technical Specification :**

Custom Gasoline Direct Injection Engine Bench



**Laboratory Instrument India**