



Email : sales@laboratoryinstrumentindia.com

Product Name :
Custom Gasoline Direct Injection Engine Bench

Product Code :
LBNY-0005-121000231



Description :

This custom trainer can be designed to incorporate any vehicle system you require. A Custom Gasoline Engine Bench equipped with an engine from any vehicle of your choice.

Standard Equipment and Features:

- All OEM components, wiring, connectors and wire colors
- Integrated ECM breakout box
- Emergency stop button
- Manual fault box for engine performance
- No component simulation
- OBDII data link connector.

Educational Advantages:

- Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation
- Demonstration of major vehicle systems to groups of students without the access limitations of a complete vehicle

Technical Specification :

Custom Gasoline Direct Injection Engine Bench



Laboratory Instrument India