



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway TXV H Block

**Product Code :**  
LBNY-0005-121000233



**Description :**

**Educational Advantage:**

Demonstrates the internal components of the TXV (H-Block) refrigerant flow control device.

The Cutaway TXV H Block cutaway trainer is designed to demonstrate the internal operation of a typical TXV H-block used to control refrigerant flow in today's HVAC systems.

**Technical Specification :**

Cutaway TXV H Block



**Laboratory Instrument India**