



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Wiper And Washer System Trainer

**Product Code :**  
LBNY-0005-121000234



**Description :**

Terminal block receptacles are provided for diagnostic test points. The Wiper And Washer System Trainer allows for the demonstration of a complete wiper and washer system using real-world components. It includes the following:

- Wiper/washer switches
- Lockable fault box with 4 faults
- Section of wiper transmission
- OEM Wiper motor
- Washer fluid pump
- Linkage enclosed in a clear plastic enclosure for safety
- Requires a 12VDC power supply to operate
- Test points for diagnosis.

**Educational Advantages:**

- Ability to test and verify a basic wiper and washer circuit.
- Four operational faults can be inserted to promote student diagnostic troubleshooting competencies.
- Ability to visualize the operation of wiper/washer switches, wiper motor and washer pump.
- Includes wiring schematics.

**Technical Specification :**

Wiper And Washer System Trainer



**Laboratory Instrument India**