



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Variable Displacement Air Conditioner
Compressor GM

Product Code :
LBNY-0005-121000235



www.laboratoryinstrumentindia.com

Description :

Features a control that allows engagement or disengagement of the compressor clutch circuit. The cutaway swash plate compressor is designed to allow students to see the internal operation of a variable displacement AC compressor. Swash plate can be manually moved to demonstrate different output modes.

Main Features:

Features a control that allows engagement or disengagement of the compressor clutch circuit
Allows visualization of internal operation of a swash plate/piston type variable displacement A/C compressor
Requires a separate power supply or 12VDC battery to operate
Driven by hand crank.

Technical Specification :

Cutaway Variable Displacement Air Conditioner Compressor GM



Laboratory Instrument India

