




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Power Window System Trainer GM Four Door

**Product Code :**  
LBNY-0005-121000236



**Description :**

Terminal block receptacles are provided for diagnostic test points. The Four-door Power Window System Trainer allows for the demonstration of a complete power window system using real world components.

**Educational Advantages:**

Electrical test points are included on the trainer for circuit diagnosis.

User manual included.

Trainer requires an external 12VDC power supply to operate.

OEM components (window motors, switches, with OEM matching wire colors).

A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.

**Technical Specification :**

Power Window System Trainer GM Four Door



**Laboratory Instrument India**