



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Variable Displacement Air Conditioning
Compressor With Slipping Clutch GM

Product Code :
LBNY-0005-121000239



www.laboratoryinstrumentindia.com

Description :

The Cutaway Variable Displacement wash plate type design is equipped with an operational 12VDC electromagnetic clutch that can be applied or disengaged with a switch. Also includes a feature that allows the demonstration of improper performance caused from a slipping compressor.

Main Features:

- Can demonstrate the improper performance caused from a slipping clutch.
- Allows visualization of variable displacement wash plate type design is equipped with an operational OEM 12VDC electromagnetic clutch.
- Driven by hand crank.
- Clutch control is protected with a 5A replaceable fuse.
- Requires a separate power supply or 12VDC battery to operate.

Technical Specification :

Cutaway Variable Displacement Air Conditioning Compressor With Slipping Clutch GM



Laboratory Instrument India

