



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Milling Machine

Product Code :
LBNY-0004-730003



Description :

Workshop Milling Machine

Technical Specification :

Milling Machine

Milling machine is defined as the type of machine in which the metal is removed by causing the job to be fed against a revolving cutting tool called the miller cutter. There are rotary cutters which covers a wide variety of different operations.

It is capable to perform its operations at small scale to larger and also heavy duty milling tasks. Now days it is commonly used for various processes in the industry and also in machine shops to provide accurate sizes and shapes of machining parts. Milling can also be done with wide range of tools.

Type: universal

Machining Capacity: Light Duty

Dimension(L*W*H): 1.4X1.6X2.2M

Weight (KG): 1000

Automatic Grade: Manual

Table Travel (X) (mm): 400 mm

Table Travel (Z) (mm): 0 mm

Range of Spindle Speed (r.p.m): 115 - 1750 r.p.m

Processing Types: Metal

Workpiece and Weight: Instrument Machine Tool

CNC or Not: Normal

Table Travel (Y) (mm): 230 mm

Workshop Milling Machine, Workshop Milling Machine Equipments, Workshop Milling Machine Tools, Workshop Milling Machine Tool Kits, Workshop Milling Machine Manufacturers, Workshop Milling Machine Suppliers from India, China, Kenya



Laboratory Instrument India