



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blower Motor System Trainer

Product Code :
LBNY-0005-121000240



Description :

Terminal block receptacles are provided for diagnostic test points. The trainer requires a 12V power supply to operate. The Blower Motor System Trainer allows for the demonstration of a complete blower motor system using real-world components.

Educational Advantages:

Four operational faults can be inserted to promote student diagnostic troubleshooting competencies.
Ability to test and verify a basic blower motor circuit.
Ability to visualize the operation of a blower motor control unit.
Includes wiring schematics.

Included Components:

Blower motor resistor.
Fault box with 4 faults.
HVAC control module.
Blower motor.
Test points for diagnosis.

Technical Specification :

Blower Motor System Trainer



Laboratory Instrument India