



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Door Lock System Trainer Four Door

Product Code :
LBNY-0005-121000241



Description :

The Power Door Lock System Trainer Four Door Hyundai is a functional Power Door Lock System Trainer mounted on a plastic training board with the following features:

- OEM components (door lock actuators, switches, with OEM matching wire colors).
- Electrical test points are included on the trainer for circuit diagnosis.
- Students can wire this circuit using supplied leads and included OEM wiring diagrams.
- A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.
- User manual included.
- Trainer requires an external 12VDC power supply to operate.

Educational Advantages:

- Fault box with 4 faults.
- Test points for diagnosis.
- Power door lock actuators (4ea).
- Door lock relay.
- Set of 21 jumper leads, 24(60cm) (6 red, 1 black, 14 blue).

Technical Specification :

Power Door Lock System Trainer Four Door



Laboratory Instrument India