



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Horn System Trainer

**Product Code :**  
LBNY-0005-121000246



**Description :**

The Horn System Trainer allows for the demonstration of a complete horn system using real-world components. The trainer requires a 12V power supply to operate. Terminal block receptacles are provided for diagnostic test points.

**Educational Advantages:**

Four operational faults can be inserted to promote student diagnostic troubleshooting competencies.  
Ability to visualize the operation of horn relays and switches.  
Ability to test and verify a basic horn circuit.  
Includes wiring schematics.

**Included Components:**

Horn  
Horn switch  
Horn relay  
Test points for diagnosis  
Fault box with 4 faults.

**Technical Specification :**

Horn System Trainer



**Laboratory Instrument India**