



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Window System Trainer GM Two Door

Product Code :
LBNY-0005-121000251



Description :

Terminal block receptacles are provided for diagnostic test points. The Two-door Power Window System Trainer allows for the demonstration of a complete power window system using real-world components.

Educational Advantages:

Ability to test and verify a basic 2-door power window circuit.

Four operational faults can be inserted to promote student diagnostic troubleshooting competencies.

Includes wiring schematics.

Ability to visualize the operation of power window actuators and switches.

Main Components:

Application: 2007+ GM Chevrolet Cobalt

Fault box with 4 faults.

Test points for diagnosis.

Power window motors (2ea).

Master and slave switches.

Dimensions: 22 x 22 x 8 in (55.9 x 55.9 x 20.3 cm)

Weight: 13 lb (5.91 kg).

Technical Specification :

Power Window System Trainer GM Two Door



Laboratory Instrument India