



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Door Lock System Trainer Two Door

Product Code :
LBNY-0005-121000252



Description :

The Power Door Lock System Trainer Two Door Hyundai is a functional Power Door Lock System Trainer mounted on a plastic training board with the following features:
A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.
Electrical test points are included on the trainer for circuit diagnosis.
User manual included.
Trainer requires an external 12VDC power supply to operate.
OEM components (door lock actuators, switches, with OEM matching wire colors).
Students can wire this circuit using supplied leads and included OEM wiring diagrams.

Advantages:

Power door lock actuators (2ea).
Test points for diagnosis.
Door lock relay.
Set of 21 jumper leads, 24 (60cm) (6 red, 1 black, 14 blue).
Fault box with 4 faults.

Technical Specification :

Power Door Lock System Trainer Two Door

Laboratory Instrument India

