



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Window System Trainer Four Door

Product Code :
LBNY-0005-121000254



www.laboratoryinstrumentindia.com

Description :

The Power Window System Trainer Four Door Hyundai mounted on a plastic training board with the following features :

Electrical test points are included on the trainer for circuit diagnosis.

Trainer requires an external 12VDC power supply to operate.

Students can wire this circuit using supplied leads and included OEM wiring diagrams.

OEM components (window motors, switches, with OEM matching wire colors).

A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.

User manual included.

Educational Advantages:

Fault box with 4 faults.

Test points for diagnosis.

Power window motors.

Master and slave switches.

Set of 21 jumper leads, 24 (60cm) (6 red, 1 black, 14 blue).

Technical Specification :

Power Window System Trainer Four Door

Laboratory Instrument India

