




Email : sales@laboratoryinstrumentindia.com

Product Name :
Coil On Plug Ignition System Trainer

Product Code :
LBNY-0005-121000255



Description :

Terminal block receptacles are provided for diagnostic test points. A 12V source is required to operate the trainer. The Coil On Plug Ignition System Trainer allows for the demonstration of a complete ignition system with distributor using real-world components.

Educational Advantages:

Ability to test and verify a basic ignition circuit without distributor
Ability to visualize the operation of a direct ignition system using coil-on-plug ignition technology
Includes wiring schematics.
Four operational faults can be inserted to promote student diagnostic troubleshooting competencies.

Included Components:

Electric motor driven crankshaft and/or camshaft sprockets
OBDII data link connector
Engine Control Module (ECM)
Crankshaft position (CKP) sensor and camshaft position (CMP) sensor
Fault box with 4 faults
Test points for diagnosis
Ignition coils
Spark plugs

Technical Specification :

Coil On Plug Ignition System Trainer



Laboratory Instrument India