



Email : sales@laboratoryinstrumentindia.com

Product Name :
Audio System Trainer

Product Code :
LBNY-0005-121000257



Description :

The Audio System Trainer is mounted on a plastic training board with OEM components. It contains the following features:

- A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.
- Electrical test points allowing student circuit diagnosis.
- Audio control unit.
- Vehicle antenna.
- Two speakers.

Educational Advantages:

- Faults can be inserted into the circuit that create real-world symptoms and abnormal operation.
- OEM components and color matching wiring are used.
- Students must connect all components using included jumper wires to reinforce relationships between components and OEM wiring diagrams.

Technical Specification :

Audio System Trainer

Laboratory Instrument India

