



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cutaway Differential With Adjustment Capabilities

**Product Code :**  
LBNY-0005-121000267



**Description :**

The Cutaway Differential With Adjustment Capabilities trainer is an application that can be chosen to demonstrate the operation of a differential equipped with adjustment screws.

**Main Features:**

Screws easily adjust the pinion depth and backlash for quick pattern demonstrations.  
Allows visualization of the internal function of each component.  
Driven by a hand crank.

**Technical Specification :**

Cutaway Differential With Adjustment Capabilities



**Laboratory Instrument India**