



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Differential With Adjustment Capabilities

Product Code :
LBNY-0005-121000267



Description :

The Cutaway Differential With Adjustment Capabilities trainer is an application that can be chosen to demonstrate the operation of a differential equipped with adjustment screws.

Main Features:

Screws easily adjust the pinion depth and backlash for quick pattern demonstrations.
Allows visualization of the internal function of each component.
Driven by a hand crank.

Technical Specification :

Cutaway Differential With Adjustment Capabilities



Laboratory Instrument India