



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Shaping Machine

**Product Code :**  
LBNY-0004-730006



**Description :**

Workshop Shaping Machine

**Technical Specification :**

Shaping Machine

The shaping machine is a plain but valuable machine. The purpose behind using this machine is to remove material such as steel or aluminum, to produce a flat surface. Moreover, it is also used to manufacture gears such as rack and pinion systems and other complex shapes.

Core Components: Gearbox, Gear

Dimension(L\*W\*H): 2000\*1300\*1500

Power (kW): 3

Name: Metal Shaping Machine

Worktable Size: 630

Weight: 1800

Shaping length:630mm

Workshop Shaping Machine, Workshop Shaping Machine Equipments, Workshop Shaping Machine Tools, Workshop Shaping Machine Tool Kits, Workshop Shaping Machine Manufacturers, Workshop Shaping Machine Suppliers from India, China, Kenya

---

# Laboratory Instrument India

