



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Shaping Machine

Product Code :
LBNY-0004-730006



Description :

Workshop Shaping Machine

Technical Specification :

Shaping Machine

The shaping machine is a plain but valuable machine. The purpose behind using this machine is to remove material such as steel or aluminum, to produce a flat surface. Moreover, it is also used to manufacture gears such as rack and pinion systems and other complex shapes.

Core Components: Gearbox, Gear

Dimension(L*W*H): 2000*1300*1500

Power (kW): 3

Name: Metal Shaping Machine

Worktable Size: 630

Weight: 1800

Shaping length:630mm

Workshop Shaping Machine, Workshop Shaping Machine Equipments, Workshop Shaping Machine Tools, Workshop Shaping Machine Tool Kits, Workshop Shaping Machine Manufacturers, Workshop Shaping Machine Suppliers from India, China, Kenya

Laboratory Instrument India

