



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
SRS Air Bag Trainer

**Product Code :**  
LBNY-0005-121000272



#### **Description :**

The SRS Air Bag Trainer allows students to understand the fundamentals of air bag systems as used in today's vehicles. Using the trainer, students will learn to recognize the main components, study various operation settings, and experiment with analytic and diagnostic methods using a diagnostic scan tool. Built with OEM components, the trainer functions in the same way as in the original vehicle with the difference that a light is used instead of the exploding airbag. A built-in lockable manual fault box is included.

#### **Educational Advantages:**

Ability to demonstrate complete SRS (Supplemental Restraint System) function in a classroom  
Steering wheel can be safely removed and reinstalled during diagnosis or system testing  
Full OEM diagnostic functionality supporting OBDII scan tool readout with OEM DTCs  
Instrumentation demonstrating air bag deployment.

#### **Equipment And Features:**

Integrated SRS breakout box with test points for all terminals sized for standard meter leads  
All OEM components, wiring, connectors and wire colors  
OBDII data link connector  
Manual fault box for SRS control module with 10 faults.

#### **Technical Specification :**

SRS Air Bag Trainer



**Laboratory Instrument India**