



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Automatic Transaxle

Product Code :
LBNY-0005-121000273



www.laboratoryinstrumentindia.com

Description :

The cutaway trainer demonstrates the internal operation of a automatic transaxle and permits the student to visualize and safely operate the system. This allows the cutaway transaxle to be easily shifted through all speeds. For training purposes the cutaway uses manually controlled air pressure to apply and release the clutches and bands.

Educational Advantages:

Gear changes can be made using regulated air pressure and provided controls (requires shop air pressure)
Driven by hand crank
Allows visualization of gear changes, planetary gear sets, clutches and bands

Technical Specification :

Dimensions: 40 x 32 x 38 in (101.5 x 81.2 x 105.4 cm)
Weight: 350 lb (159 kg).



Laboratory Instrument India

