



Email : sales@laboratoryinstrumentindia.com

Product Name :
Engine Stand Starting Station

Product Code :
LBNY-0005-121000276



Description :

Advantages:

Provides a single starting station that can be used with multiple engine stands
All rotating components are covered with safety shields or guards to prevent injuries
Accessibility of engine components is enhanced by not having the requirement of whole vehicles in the shop
Continuity of subject matter is provided by students working on identical engines
Accessibility of engine components is enhanced by not having the requirement of whole vehicles in the shop
Ability to verify proper engine systems operation before and after engine repair.

Applications:

Custom design of engine stand to meet individual requirements
Gasoline or diesel engine of your choice

Standard Equipment And Features:

Includes the custom engine stand and the custom starting station available separately
Integrated PCM breakout box with test points for all terminals sized for standard meter leads
Custom machined mounting plate for engine mounting to stand
All OEM components, wiring, connectors and wire colors
No component simulation
OBDII data link connector.

Technical Specification :

Engine Stand Starting Station



Laboratory Instrument India