



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Engine Stand Starting Station

**Product Code :**  
LBNY-0005-121000276



#### **Description :**

##### **Advantages:**

Provides a single starting station that can be used with multiple engine stands  
All rotating components are covered with safety shields or guards to prevent injuries  
Accessibility of engine components is enhanced by not having the requirement of whole vehicles in the shop  
Continuity of subject matter is provided by students working on identical engines  
Accessibility of engine components is enhanced by not having the requirement of whole vehicles in the shop  
Ability to verify proper engine systems operation before and after engine repair.

##### **Applications:**

Custom design of engine stand to meet individual requirements  
Gasoline or diesel engine of your choice

##### **Standard Equipment And Features:**

Includes the custom engine stand and the custom starting station available separately  
Integrated PCM breakout box with test points for all terminals sized for standard meter leads  
Custom machined mounting plate for engine mounting to stand  
All OEM components, wiring, connectors and wire colors  
No component simulation  
OBDII data link connector.

#### **Technical Specification :**

Engine Stand Starting Station



**Laboratory Instrument India**