



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway Supercharger

**Product Code :**  
LBNY-0005-121000277



**Description :**

A Cutaway Supercharger will allow your students to rotate the supercharger by hand to visualize the operating principals that make them function.

**Advantages:**

Provides students with full accessibility to observe all internal components of a supercharger and how they interact with each other  
Built to application of your choice.  
Able to be rotated by hand for demonstration purposes.

**Technical Specification :**

Cutaway Supercharger



**Laboratory Instrument India**