



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Supercharger

Product Code :
LBNY-0005-121000277



Description :

A Cutaway Supercharger will allow your students to rotate the supercharger by hand to visualize the operating principals that make them function.

Advantages:

Provides students with full accessibility to observe all internal components of a supercharger and how they interact with each other
Built to application of your choice.
Able to be rotated by hand for demonstration purposes.

Technical Specification :

Cutaway Supercharger



Laboratory Instrument India