



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Timing System Belt Trainer

**Product Code :**  
LBNY-0005-121000278



**Description :**

The Timing System Belt Trainer visually identifies all of the components in a typical belt type camshaft drive system for OHC engines.

**Advantages:**

Allows full access of all OHC drive belt components without limitations of vehicle components, covers and shrouds blocking view  
Students are able to have clear access to all components for skill practice prior to working on a complete vehicle.  
Can be used for clear demonstrations of common service inspections and repair procedures of:

- a. Drive belt
- b. Cam and Crank oil seals
- c. Cam and Crank gears
- d. Hydraulic belt tensioner
- e. Water pump drive gear
- f. Idler puller
- g. Tensioner pulley.

**Technical Specification :**

Timing System Belt Trainer

---

# Laboratory Instrument India

