



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Workshop Movable Crane

**Product Code :**  
LBNY-0004-730007



**Description :**

Workshop Movable Crane

**Technical Specification :**

Movable Crane  
Rated Lifting Moment: 2000kg  
Application: Workshop,  
Rated Lifting Moment: 2000kg  
Span: 3000mm  
Control Method: Hand push or electric controlled  
Crane type: Single Beam Light Duty Gantry Crane  
Crane feature: Adjustable height  
Rated Loading Capacity: 3000kg  
Max. Lifting Height: 4000mm  
Material: Structural Steel or Aluminum Alloy  
Power Source: 3 Phase 380V 50hz  
Lifting mechanism: Electric Chain Hoist

Workshop Movable Crane, Workshop Movable Crane Equipments, Workshop Movable Crane Tools, Workshop Movable Crane Tool Kits, Workshop Movable Crane Manufacturers, Workshop Movable Crane Suppliers from India, China, Kenya



**Laboratory Instrument India**