



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Movable Crane

Product Code :
LBNY-0004-730007



Description :

Workshop Movable Crane

Technical Specification :

Movable Crane
Rated Lifting Moment: 2000kg
Application: Workshop,
Rated Lifting Moment: 2000kg
Span: 3000mm
Control Method: Hand push or electric controlled
Crane type: Single Beam Light Duty Gantry Crane
Crane feature: Adjustable height
Rated Loading Capacity: 3000kg
Max. Lifting Height: 4000mm
Material: Structural Steel or Aluminum Alloy
Power Source: 3 Phase 380V 50hz
Lifting mechanism: Electric Chain Hoist

Workshop Movable Crane, Workshop Movable Crane Equipments, Workshop Movable Crane Tools, Workshop Movable Crane Tool Kits, Workshop Movable Crane Manufacturers, Workshop Movable Crane Suppliers from India, China, Kenya



Laboratory Instrument India