



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Serpentine System Trainer

**Product Code :**  
LBNY-0005-121000280



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

This configuration is very common in modern vehicle belt drive applications. The Serpentine System Belt Trainer is a visual representation of an accessory drive system using both sides of the drive belt to transmit torque.

**Advantages:**

Allows student practice for replacement and normal service procedures of:

- a. Tensioner
- b. Drive belt
- c. Idler pulley

Provides accessibility for teaching the proper methods of belt and tensioner condition inspection without limitations of having a complete vehicle with shrouds, components and covers blocking access.

An industry standard belt wear gauge is included.

**Technical Specification :**

Serpentine System Trainer



**Laboratory Instrument India**

