



Email : sales@laboratoryinstrumentindia.com

Product Name :
2.2l Engine Bench With HVAC

Product Code :
LBNY-0005-121000281



www.laboratoryinstrumentindia.com

Description :

The trainer is easy to work on and all major components, wiring and connectors match OEM configuration. 2.2l Engine Bench With HVAC Trainer is mounted on a mobile stand constructed of 2 (5cm²) tubular steel equipped with four 5 phenol casters (2-fixed and 2-swivel w/locks). The 2.2l Engine Bench With HVAC is a fully operational engine mounted on a mobile stand. Trainer is based on 2007 and newer 2.2L global engine platform. Trainer is equipped with safety guards covering all rotating components.

Advantages:

Engine and A/C systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation
Visually demonstrates refrigerant change of state
Hands on demonstration of electric assist power steering
Demonstration of systems to groups of students without the access limitations of a complete vehicle

Standard Equipment and Features:

Safety covers for all rotating components
Fuel tank
12V battery included with automatic charger
Automatic transaxle
Refrigerant high and low-side pressure gauges
Exhaust system w/ catalytic converter
Manual fault box for engine performance
HD square tubing mobile stand with phenol casters
Emergency stop button.

Technical Specification :

2.2l Engine Bench With HVAC



Laboratory Instrument India