



Email : sales@laboratoryinstrumentindia.com

Product Name :
Swivel Engine Stand No Engine

Product Code :
LBNY-0005-121000282



www.laboratoryinstrumentindia.com

Description :

Custom machined mounting plate. Sturdy Swivel Stand used for mounting engine to meet your specific needs in engine repair.

Features:

Sturdy construction
Swivel handle
Custom machines mounting plate for engine mounting to stand.

Technical Specification :

Swivel Engine Stand No Engine



Laboratory Instrument India