



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Swivel Engine Stand No Engine

**Product Code :**  
LBNY-0005-121000282



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Custom machined mounting plate. Sturdy Swivel Stand used for mounting engine to meet your specific needs in engine repair.

**Features:**

Sturdy construction  
Swivel handle  
Custom machines mounting plate for engine mounting to stand.

**Technical Specification :**

Swivel Engine Stand No Engine



**Laboratory Instrument India**