



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
1.6l Engine Bench

**Product Code :**  
LBNY-0005-121000284



#### **Description :**

Trainer is mounted on a mobile stand constructed of 2(5cm) tubular steel equipped with four 5(12.5cm) phenol casters (2-fixed and 2-swivel w/locks). 1.6l Engine Bench Trainer is equipped with safety guards covering all rotating components. The engine performance trainer is a fully operational engine mounted on a mobile stand. The trainer is easy to work on and all major components, wiring and connectors match OEM configuration. Trainer is based on 2007 and newer 1.6L global engine platform.

#### **Advantages:**

Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation.

Demonstration of systems to groups of students without the access limitations of a complete vehicle.

#### **Features:**

Sequential fuel injection w/coil-on-plug

1.6L DOHC engine

Vehicle OEM VIN number supplied with trainer

Adaptive charging system

#### **Technical Specification :**

Fuel requirements: 87 octane gasoline

Dimensions: 29 W x 48 L x 50 H (72.5 x 120 x 125 cm) w/o pallet 35 W x 60 L x 56 H (90 x 150 x 140 cm) w/pallet

Weight: 775 lbs (340.9 kg) w/o pallet 850 lbs (386.4 kg) w/pallet.



**Laboratory Instrument India**